

**Box and Trays**  
**Caja y Estuches**

Ref.	Box and Trays   Caja y Estuches
09.50.02.00180	Box, 250 x 435 x 130mm, DHS-DCS 4.5/5.0mm <i>Caja, DHS-DCS 4.5/5.0mm</i>
09.50.06.00180	Instrument Tray, Set-1, DHS-DCS 4.5/5.0mm <i>Bandeja Instrumentos, Set-1, DHS-DCS 4.5/5.0mm</i>
09.50.06.00181	Instrument Tray, Set-2, DHS-DCS 4.5/5.0mm <i>Bandeja Instrumentos, Set-2, DHS-DCS 4.5/5.0mm</i>
09.50.06.00182	Implant Tray, Set-1, DHS-DCS 4.5/5.0mm <i>Bandeja Implantes, Set-1, DHS-DCS 4.5/5.0mm</i>
09.50.06.00183	Implant Tray, Set-2, DHS-DCS 4.5/5.0mm <i>Bandeja Implantes, Set-2, DHS-DCS 4.5/5.0mm</i>
09.50.10.00150	Dispenser for K-Wire 150mm <i>Dispensador para Agujas de Kirschner 150mm</i>
09.50.10.00230	Dispenser for K-Wire 230mm <i>Dispensador para Agujas de Kirschner 230mm</i>



**PLATES | PLACAS**

**DCS Locking Plate, 95 Degrees**  
**Placa Bloqueo, DCS, 95 Grados**

Ref.	Nº Holes   Nº Agujeros	 
27.02.00.09506	06	
27.02.00.09508	08	
27.02.00.09510	10	
27.02.00.09512	12	
27.02.00.09514	14	



**DHS Locking Plate, 135 Degrees, Std Barrel**  
**Placa Bloqueo DHS, 135 Grados, Barril Std**

Ref.	Nº Holes   Nº Agujeros	 
27.04.00.13502	02	
27.04.00.13503	03	
27.04.00.13504	04	
27.04.00.13505	05	
27.04.00.13506	06	
27.04.00.13508	08	
27.04.00.13510	10	
27.04.00.13512	12	
27.04.00.13514	14	



## 4.5mm/5.0mm DHS | DCS

### DHS Locking Plate, 135 Degrees, Short Barrel *Placa Bloqueo DHS, 135 Grados, Barril Corto*

Ref.	Nº Holes   Nº Agujeros	 
27.04.01.13502	02	
27.04.01.13503	03	
27.04.01.13504	04	
27.04.01.13505	05	
27.04.01.13506	06	




### Ø12,7 DHS-DCS Lag Screw, 22mm Threaded *Tornillo DHS-DCS, Rosca 22mm*

Ref.	Length   Longitud
15.22.22.12050	50mm
15.22.22.12055	55mm
15.22.22.12060	60mm
15.22.22.12065	65mm
15.22.22.12070	70mm
15.22.22.12075	75mm
15.22.22.12080	80mm
15.22.22.12085	85mm
15.22.22.12090	90mm
15.22.22.12095	95mm
15.22.22.12100	100mm
15.22.22.12105	105mm
15.22.22.12110	110mm
15.22.22.12115	115mm
15.22.22.12120	120mm
15.22.22.12125	125mm
15.22.22.12130	130mm



### DHS-DCS Compression Screw *Tornillo de Compresión DHS-DCS*

Ref.	Length   Longitud	
15.24.01.00000	DHS-DCS Compression Screw, Ti	



# 4.5mm/5.0mm DHS | DCS



## INSTRUMENTS INSTRUMENTOS

09.01.01.20100	Drill Bit, Ø2.0 x 100mm, twist length 30mm, AO Coupling <i>Broca Ø2.0 x 100mm, H-30, AO</i>	2
09.01.01.45195	Drill Bit, Ø4.5 x 195mm, twist length 50mm, AO Coupling <i>Broca Ø4.5 x 195mm, H-50, AO</i>	2
09.01.14.32195	Drill Bit, Ø3.2 x 195mm, twist length 50mm, AO Coupling, Yellow Code <i>Broca Ø3.2 x 195mm, H-50, AO, Ama</i>	2
09.01.14.43195	Drill Bit, Ø4.3 x 195mm, twist length 50mm, AO Coupling, Green Code <i>Broca Ø4.3 x 195mm, H-50, AO, Gri</i>	2
09.03.01.45000	Tap, Cortez, Ø4.5 x 1.75mm, AO Coupling <i>Macho Cortical Ø4.5 x 1.75mm, AO</i>	1
09.03.02.65000	Tap, Cancellous, Ø6.5 x 2.75mm, AO Coupling <i>Macho Esponjosa Ø6.5 x 2.75mm, AO</i>	1
09.03.05.12500	Tap, DHS-DCS, Ø12.5mm <i>Macho DHS-DCS Ø12.5mm</i>	1
09.04.02.00000	T-Handle, DHS-DCS <i>Mango en T, DHS-DCS</i>	1
09.04.04.15020	Handle, Cannulated, AO Coupling, 150mm, Yellow <i>Mango Canulado, AO, 150mm, Ama</i>	1
09.04.07.00020	T-Handle, Cannulated, AO Coupling, Yellow <i>Mango en T Canulado, AO, Ama</i>	1
09.05.00.45000	Drill Guide, Ø4.5/Ø3.2mm, Universal <i>Guía de Perforación Ø4.5/Ø3.2mm, Universal</i>	1
09.05.02.45320	Drill Guide, Ø4.5/Ø3.2mm, Double <i>Guía de Perforación Ø4.5/Ø3.2mm, Doble</i>	1
09.05.04.00095	Drill Guide, DCS - 95° <i>Guía de Perforación DCS - 95°</i>	1
09.05.04.00135	Drill Guide, DHS - 135° <i>Guía de Perforación DHS - 135°</i>	1
09.05.05.43000	Drill Guide, Ø4.3mm, Threaded <i>Guía de Perforación Ø4.3mm, Roscado</i>	2
09.06.00.32000	Drill Sleeve, Ø3.2mm <i>Manga de Perforación Ø3.2mm</i>	1
09.06.04.45000	Holding Sleeve, Ø4.5mm and Larger <i>Manga de Sujeción Ø4.5mm y Superiores</i>	1
09.06.05.00000	Centering Sleeve, DHS-DCS <i>Manga Centradora, DHS-DCS</i>	1
09.07.05.25170	Shaft Screwdriver, Torx-25, 170mm, AO Coupling, Barrel Ø7.0mm <i>Destornillador, TX-25, 170mm, AO, C-7.0</i>	2
09.08.01.00120	Depth Gauge, 120mm <i>Medidor de Profundidad 120mm</i>	1
09.08.02.25160	Gauge 160mm, Direct Measuring, for Wire 230mm, Ø up to 2.5mm <i>Medidor 90mm, Aguja Ø2.0 x 230mm</i>	1
09.18.01.00040	Torque Limiting Attachment, 4.0 Nm, AO Coupling <i>Acoplamiento/Limitador de Torque, 4.0 Nm, AO</i>	1
09.19.01.00020	Tensioner, Span 20mm, Articulated <i>Tensor, 20mm, Articulado</i>	1
09.20.01.11000	Socket Wrench 11mm, Universal Joint <i>Llave de Tubo 11mm, Junta Universal</i>	1
09.20.02.00000	DHS-DCS Insertion Wrench <i>Llave de Inserción, DHS-DCS</i>	1
09.21.03.00001	Reamer, DCS Complete <i>Fresa Completa, DCS</i>	1
09.21.03.00002	Reamer, DHS Complete, Short <i>Fresa Completa, DHS, Corto</i>	1
09.21.03.00003	Reamer, DHS Complete, Long <i>Fresa Completa, DHS, Largo</i>	1
09.22.02.00010	Coupling Screw, DHS-DCS, Cannulated <i>Tornillo de Acoplamiento, DHS-DCS, Canulado</i>	1
09.23.01.00000	Impactor, DHS-DCS <i>Impactor, DHS-DCS</i>	1
09.41.01.00025	Screw Extractor, for Torx-25, AO Coupling <i>Extractor de Tornillo, TX-25, AO</i>	1

## DHS-DCS LAG SCREWS TORNILLOS DHS-DCS

15.22.22.12050	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 50mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 50mm</i>	1
15.22.22.12055	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 55mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 55mm</i>	1
15.22.22.12060	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 60mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 60mm</i>	1
15.22.22.12065	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 65mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 65mm</i>	1
15.22.22.12070	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 70mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 70mm</i>	1
15.22.22.12075	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 75mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 75mm</i>	1
15.22.22.12080	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 80mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 80mm</i>	1
15.22.22.12085	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 85mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 85mm</i>	1

## DHS-DCS LAG SCREWS TORNILLOS DHS-DCS

15.22.22.12090	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 90mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 90mm</i>	1
15.22.22.12095	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 95mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 95mm</i>	1
15.22.22.12100	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 100mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 100mm</i>	1
15.22.22.12105	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 105mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 105mm</i>	1
15.22.22.12110	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 110mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 110mm</i>	1
15.22.22.12115	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 115mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 115mm</i>	1
15.22.22.12120	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 120mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 120mm</i>	1
15.22.22.12125	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 125mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 125mm</i>	1
15.22.22.12130	DHS-DCS Lag Screw, Ti, 22mm Threaded, Ø12.7 x 130mm <i>Tornillo DHS-DCS, Ti, Rosca 22mm, Ø12.7 x 130mm</i>	1

## DHS-DCS COMPRESSION SCREW TORNILLO DE COMPRESIÓN DHS-DCS

15.24.01.00000	DHS-DCS Compression Screw, Ti <i>Tornillo de Compresión DHS-DCS, Ti</i>	5
----------------	--	---

## DCS LOCKING PLATE PLACA BLOQUEO DCS

27.02.00.09506	DCS Locking Plate, Ti, 95 Degrees, 06 Holes <i>Placa Bloqueo, Ti, DCS, 95 Grados, 6-A</i>	2
27.02.00.09508	DCS Locking Plate, Ti, 95 Degrees, 08 Holes <i>Placa Bloqueo, Ti, DCS, 95 Grados, 8-A</i>	2
27.02.00.09510	DCS Locking Plate, Ti, 95 Degrees, 10 Holes <i>Placa Bloqueo, Ti, DCS, 95 Grados, 10-A</i>	2
27.02.00.09512	DCS Locking Plate, Ti, 95 Degrees, 12 Holes <i>Placa Bloqueo, Ti, DCS, 95 Grados, 12-A</i>	1
27.02.00.09514	DCS Locking Plate, Ti, 95 Degrees, 14 Holes <i>Placa Bloqueo, Ti, DCS, 95 Grados, 14-A</i>	1

## DHS LOCKING PLATE PLACA BLOQUEO DHS

27.04.00.13502	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 02 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 2-A</i>	2
27.04.00.13503	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 03 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 3-A</i>	2
27.04.00.13504	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 04 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 4-A</i>	2
27.04.00.13505	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 05 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 5-A</i>	2
27.04.00.13506	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 06 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 6-A</i>	2
27.04.00.13508	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 08 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 8-A</i>	2
27.04.00.13510	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 10 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 10-A</i>	1
27.04.00.13512	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 12 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 12-A</i>	1
27.04.00.13514	DHS Locking Plate, Ti, 135 Degrees, Std Barrel, 14 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Std, 14-A</i>	1
27.04.01.13502	DHS Locking Plate, Ti, 135 Degrees, Short Barrel, 02 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Corto, 2-A</i>	1
27.04.01.13503	DHS Locking Plate, Ti, 135 Degrees, Short Barrel, 03 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Corto, 3-A</i>	1
27.04.01.13504	DHS Locking Plate, Ti, 135 Degrees, Short Barrel, 04 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Corto, 4-A</i>	1
27.04.01.13505	DHS Locking Plate, Ti, 135 Degrees, Short Barrel, 05 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Corto, 5-A</i>	1
27.04.01.13506	DHS Locking Plate, Ti, 135 Degrees, Short Barrel, 06 Holes <i>Placa Bloq DHS, Ti, 135 Grados, Barril Corto, 6-A</i>	1