



# Astrolabe

**LIFE AND MOBILITY**

## Box Components


### Humerus Proximal Full Set

Box Components Kit

09.50.03.00166


Component Trays

09.50.06.00152






Humerus Proximal  
Full Set  
3,5



# Humerus Proximal Full Set




N#	REF	DESC	QTY
	15.04.01.35012	Cortical Screw, Ti, Ø3,5 x 12mm Torx	5
	15.04.01.35014	Cortical Screw, Ti, Ø3,5 x 14mm Torx	5
	15.04.01.35016	Cortical Screw, Ti, Ø3,5 x 16mm Torx	5
	15.04.01.35018	Cortical Screw, Ti, Ø3,5 x 18mm Torx	5
	15.04.01.35020	Cortical Screw, Ti, Ø3,5 x 20mm Torx	5
	15.04.01.35022	Cortical Screw, Ti, Ø3,5 x 22mm Torx	5
	15.04.01.35024	Cortical Screw, Ti, Ø3,5 x 24mm Torx	5
	15.04.01.35026	Cortical Screw, Ti, Ø3,5 x 26mm Torx	5
	15.04.01.35028	Cortical Screw, Ti, Ø3,5 x 28mm Torx	5
	15.04.01.35030	Cortical Screw, Ti, Ø3,5 x 30mm Torx	5
	15.04.01.35032	Cortical Screw, Ti, Ø3,5 x 32mm Torx	5
	15.04.01.35034	Cortical Screw, Ti, Ø3,5 x 34mm Torx	5
	15.04.01.35036	Cortical Screw, Ti, Ø3,5 x 36mm Torx	5
1	15.04.01.35038	Cortical Screw, Ti, Ø3,5 x 38mm Torx	5
	15.04.01.35040	Cortical Screw, Ti, Ø3,5 x 40mm Torx	5
	15.04.01.35042	Cortical Screw, Ti, Ø3,5 x 42mm Torx	5
	15.04.01.35044	Cortical Screw, Ti, Ø3,5 x 44mm Torx	5
	15.04.01.35045	Cortical Screw, Ti, Ø3,5 x 45mm Torx	5
	15.04.01.35046	Cortical Screw, Ti, Ø3,5 x 46mm Torx	5
	15.04.01.35048	Cortical Screw, Ti, Ø3,5 x 48mm Torx	5
	15.04.01.35050	Cortical Screw, Ti, Ø3,5 x 50mm Torx	5
	15.04.01.35052	Cortical Screw, Ti, Ø3,5 x 52mm Torx	5
	15.04.01.35054	Cortical Screw, Ti, Ø3,5 x 54mm Torx	5
	15.04.01.35055	Cortical Screw, Ti, Ø3,5 x 55mm Torx	5
	15.04.01.35056	Cortical Screw, Ti, Ø3,5 x 56mm Torx	5
	15.04.01.35058	Cortical Screw, Ti, Ø3,5 x 58mm Torx	5
	15.04.01.35060	Cortical Screw, Ti, Ø3,5 x 60mm Torx	5






2



3

N#	REF	DESC	QTY
	15.14.01.35012	Locking Screw, Ti, Ø3,5 x 12mm Torx	5
	15.14.01.35014	Locking Screw, Ti, Ø3,5 x 14mm Torx	5
	15.14.01.35016	Locking Screw, Ti, Ø3,5 x 16mm Torx	5
	15.14.01.35018	Locking Screw, Ti, Ø3,5 x 18mm Torx	5
	15.14.01.35020	Locking Screw, Ti, Ø3,5 x 20mm Torx	5
	15.14.01.35022	Locking Screw, Ti, Ø3,5 x 22mm Torx	5
	15.14.01.35024	Locking Screw, Ti, Ø3,5 x 24mm Torx	5
	15.14.01.35026	Locking Screw, Ti, Ø3,5 x 26mm Torx	5
	15.14.01.35028	Locking Screw, Ti, Ø3,5 x 28mm Torx	5
	15.14.01.35030	Locking Screw, Ti, Ø3,5 x 30mm Torx	5
	15.14.01.35032	Locking Screw, Ti, Ø3,5 x 32mm Torx	5
	15.14.01.35034	Locking Screw, Ti, Ø3,5 x 34mm Torx	5
	15.14.01.35036	Locking Screw, Ti, Ø3,5 x 36mm Torx	5
	15.14.01.35038	Locking Screw, Ti, Ø3,5 x 38mm Torx	5
	15.14.01.35040	Locking Screw, Ti, Ø3,5 x 40mm Torx	5
	15.14.01.35042	Locking Screw, Ti, Ø3,5 x 42mm Torx	5
	15.14.01.35044	Locking Screw, Ti, Ø3,5 x 44mm Torx	5
	15.14.01.35045	Locking Screw, Ti, Ø3,5 x 45mm Torx	5
	15.14.01.35046	Locking Screw, Ti, Ø3,5 x 46mm Torx	5
	15.14.01.35048	Locking Screw, Ti, Ø3,5 x 48mm Torx	5
	15.14.01.35050	Locking Screw, Ti, Ø3,5 x 50mm Torx	5
	15.14.01.35052	Locking Screw, Ti, Ø3,5 x 52mm Torx	5
	15.14.01.35054	Locking Screw, Ti, Ø3,5 x 54mm Torx	5
	15.14.01.35055	Locking Screw, Ti, Ø3,5 x 55mm Torx	5
	15.14.01.35056	Locking Screw, Ti, Ø3,5 x 56mm Torx	5
	15.14.01.35058	Locking Screw, Ti, Ø3,5 x 58mm Torx	5
	15.14.01.35060	Locking Screw, Ti, Ø3,5 x 60mm Torx	5
	15.08.01.40012	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 12mm Torx	2
	15.08.01.40014	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 14mm Torx	2
	15.08.01.40016	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 16mm Torx	2
	15.08.01.40018	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 18mm Torx	2
	15.08.01.40020	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 20mm Torx	2
	15.08.01.40022	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 22mm Torx	2
	15.08.01.40024	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 24mm Torx	2
	15.08.01.40026	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 26mm Torx	2
	15.08.01.40028	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 28mm Torx	2
	15.08.01.40030	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 30mm Torx	2
	15.08.01.40032	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 32mm Torx	2
	15.08.01.40034	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 34mm Torx	2
	15.08.01.40036	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 36mm Torx	2
	15.08.01.40038	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 38mm Torx	2
	15.08.01.40040	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 40mm Torx	2
	15.08.01.40045	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 45mm Torx	2
	15.08.01.40050	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 50mm Torx	2
	15.08.01.40055	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 55mm Torx	2
	15.08.01.40060	Cancellous Screw, Fully Threaded, Ti, Ø4,0 x 60mm Torx	2

# Humerus Proximal Full Set













	N#	REF	DESC	QTY
	4	09.05.00.35000	Drill Guide, Ø3.5/Ø2.5mm, Universal	1
	5	09.05.00.27000	Drill Guide, Ø2.7/Ø2.0mm, Universal	1
	6	09.15.03.35008	Bending Template, for Humerus Proximal Plate, 8-H Long	1
	7	09.15.03.35005	Bending Template, for Humerus Proximal Plate, 5-H	1
	8	09.15.03.35003	Bending Template, for Humerus Proximal Plate, 3-H	1
	9	26.08.11.00003	Special Locking Plate, Ti, Humerus, Proximal, 03 Holes	1
	10	26.08.11.00005	Special Locking Plate, Ti, Humerus, Proximal, 05 Holes	1



N#	REF	DESC	QTY
11	26.08.10.00005	Special Locking Plate, Ti, Humerus, Proximal, 05 Holes - Long	1
12	26.08.10.00006	Special Locking Plate, Ti, Humerus, Proximal, 06 Holes - Long	1
13	26.08.10.00008	Special Locking Plate, Ti, Humerus, Proximal, 08 Holes - Long	1
14	26.08.10.00010	Special Locking Plate, Ti, Humerus, Proximal, 10 Holes - Long	1
15	26.08.10.00012	Special Locking Plate, Ti, Humerus, Proximal, 12 Holes - Long	1
16	09.05.06.00000	Insertion Guide, Proximal Humerus Plate	1
17	09.05.06.00010	Insertion Guide, Proximal Humerus Plate, with Nose	1











# Humerus Proximal Full Set



	N#	REF	DESC	QTY
	18	09.37.00.00000	Depth Probe, Humerus Plate	1
	19	09.06.04.27350	Holding Sleeve, Ø2.7/3.5mm	1
	20	09.17.01.06000	Periosteal Elevator, 6mm width, Curved Blade, Round Edges	1
	21	09.03.01.35000	Tap, Cortex, Ø3.5 x 1.25mm, AO Coupling	1
	22	09.03.02.40000	Tap, Cancellous, Ø4.0 x 1.75mm, AO Coupling	1
	23	09.07.05.15117	Shaft Screwdriver, Torx-15, 110mm, AO Coupling, Barrel Ø5.0mm, Blue-Red Code	2
	24	09.01.12.28125	Drill Bit, Ø2.8 x 125mm, AO Coupling, Barrel Ø6.5mm	2
	25	09.04.04.13510	Handle, Cannulated, AO Coupling, 135mm, Blue	1
	26	09.04.04.12030	Handle, Cannulated, AO Coupling, 120mm, Green	1
	27	09.01.13.28165	Graduated Drill Bit, Ø2.8 x 165mm, AO Coupling, with Drill Stop, Red Code	2

**Additional Components**

REF	Desc
09.18.01.00015	Torque Limiting Attachment, 1.5 Nm, AO

N#	REF	DESC	QTY
 28	09.01.01.25165	Drill Bit, Ø2.5 x 165mm, twist length 40mm, AO Coupling	2
 29	09.01.01.20125	Drill Bit, Ø2.0 x 125mm, twist length 30 mm, AO Coupling	2
 30	09.28.00.26160	Pulling Device Ø2.6 x 160mm, Without Nut	1
 31	09.38.02.00000	Wing Nut, for Pull Reduction Device	1
 32	09.06.09.53000	Outer Sleeve, Ø5.3mm	1
 33	09.06.00.28000	Drill Sleeve, Ø2.8mm	1
 34	09.06.02.16000	Wire Sleeve, Ø1.6mm	1
 35	09.06.01.28000	Drill Sleeve, Ø2.8mm, Threaded	2
 36	09.08.01.00060	Depth Gauge, 60mm	1
	09.50.10.00150	Dispenser for K-Wire 150mm	1
 37	17.02.03.16150	Kirschner Wire, SS, Threaded Tip, Ø1.6 x 150mm	10







AST.BOX.043  
REV.02  
02/2026