


Cannulated Bone Screws


Wide range of sizes for all your needs

AstroCAN

Astrolabe's Cannulated screw family for small and large bone fixation.



- Titanium alloy (Ti-6Al-4V);
- Various diameters and lengths;
- Self-tapping design;
- Low profile head for reduction of soft tissue damage;
- Long and short thread options;
- Torx screws driving recess;



Cannulated Screws
Ø2,7



Cannulated Screws
Ø3,5



Cannulated Screws
Ø4,0



Cannulated Screws
Ø4,5



Cannulated Screws
Ø7,0

Cannulated Conical Head Screws
Ø2,7



Specific for forefoot procedures

Cannulated screw washers per system.

Main Indications	CBS Ø2,7	CBS Ø3,5	CBS Ø4,0	CBS Ø4,5	CBS Ø7,0
Proximal Humerus		X	X		
Distal Humerus	X	X			
Proximal Radius	X				
Radial & Ulnar Styloid					
Scaphoid					
Carpal Arthrodesis					
Femoral Condyles				X	X
Patella		X	X		
Tibial Plateau			X	X	X
MTP Arthrodesis & Hallux Valgus	X				
Midfoot Injuries & Arthrodesis	X	X			
Rearfoot Injuries & Osteotomies				X	X


Cannulated Compression Screws

Wide range of sizes for all your needs

Astrolabe's Cannulated screw family for small and large bone fixation.



- Titanium alloy (Ti-6Al-4V);
- Various diameters and lengths;
- Self-tapping design;
- Headless design for reduction of soft tissue damage;
- Torx screws driving recess;



Main Indications	CCS Ø2,5	CCS Ø3,0	CCS Ø3,5	CCS Ø6,0
Proximal Humerus		X		
Distal Humerus	X	X		
Proximal Radius	X	X		
Radial & Ulnar Styloid	X	X		
Scaphoid	X	X		
Carpal Arthrodesis	X	X		
Femoral Condyles				X
Patella		X	X	
Tibial Plateau			X	X
MTP Arthrodesis & Hallux Valgus	X	X		
Midfoot Injuries & Arthrodesis	X	X	X	X
Rearfoot Injuries & Osteotomies			X	X